

# Validation Document for TETPOR AIR Pharmaceutical Grade Cartridge and DEMICAP Capsule Filters



**domnick hunter**

*PROCESS FILTRATION*

## Contents

|           |                                                                                                                                                                                   |           |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| <b>1</b>  | <b>INTRODUCTION .....</b>                                                                                                                                                         | <b>3</b>  |
| <b>2</b>  | <b>QUALITY POLICY .....</b>                                                                                                                                                       | <b>4</b>  |
| 2.1       | QUALITY ASSURANCE .....                                                                                                                                                           | 4         |
| <b>3</b>  | <b>PRODUCT TRACEABILITY.....</b>                                                                                                                                                  | <b>5</b>  |
| <b>4</b>  | <b>SUMMARY OF TEST METHODS.....</b>                                                                                                                                               | <b>6</b>  |
| <b>5</b>  | <b>PREPARATION OF PHARMACEUTICAL GRADE FILTERS .....</b>                                                                                                                          | <b>7</b>  |
| 5.1       | PHARMACEUTICAL GRADE STANDARD.....                                                                                                                                                | 7         |
| 5.2       | QUALITY OF PURIFIED WATER USED IN THE PREPARATION OF PHARMACEUTICAL GRADE FILTERS .....                                                                                           | 7         |
| <b>6</b>  | <b>CERTIFICATE OF CONFORMANCE.....</b>                                                                                                                                            | <b>8</b>  |
| <b>7</b>  | <b>PRODUCT CODING AND RANGE FOR TETPOR AIR CARTRIDGES.....</b>                                                                                                                    | <b>9</b>  |
| <b>8</b>  | <b>PRODUCT CODING AND RANGE FOR TETPOR AIR CAPSULE FILTERS.....</b>                                                                                                               | <b>9</b>  |
| <b>9</b>  | <b>ENDCAP CONFIGURATION DIAGRAMS.....</b>                                                                                                                                         | <b>10</b> |
| <b>10</b> | <b>CAPSULE DIMENSIONS.....</b>                                                                                                                                                    | <b>11</b> |
| <b>11</b> | <b>PRODUCT SPECIFICATION.....</b>                                                                                                                                                 | <b>12</b> |
| 11.1      | APPLICATION .....                                                                                                                                                                 | 12        |
| 11.2      | MATERIALS OF CONSTRUCTION FOR TETPOR AIR CARTRIDGES AND CAPSULES.....                                                                                                             | 12        |
| 11.3      | FLOW CHARACTERISTICS.....                                                                                                                                                         | 13        |
| 11.4      | OPERATING TEMPERATURES AND PRESSURES .....                                                                                                                                        | 14        |
| 11.5      | STEAM STERILISATION .....                                                                                                                                                         | 14        |
|           | <i>Autoclave.....</i>                                                                                                                                                             | <i>14</i> |
|           | <i>Steam in Place (SIP).....</i>                                                                                                                                                  | <i>14</i> |
| 11.6      | TESTS FOR BIOCOMPATIBILITY .....                                                                                                                                                  | 15        |
| 11.7      | EXTRACTABLES .....                                                                                                                                                                | 15        |
|           | <i>Test Method (1).....</i>                                                                                                                                                       | <i>15</i> |
|           | <i>Test Method (2).....</i>                                                                                                                                                       | <i>15</i> |
| 11.8      | INTEGRITY TESTING DATA .....                                                                                                                                                      | 16        |
| <b>12</b> | <b>CORRELATION OF NON-DESTRUCTIVE INTEGRITY TESTING TO LIQUID BACTERIAL CHALLENGE WITH <i>BREVUNDIMONAS DIMINUTA</i> (ATCC 19146) FOR 0.20 µM STERILISING GRADE FILTERS .....</b> | <b>16</b> |
| 12.1      | LIQUID BACTERIAL CHALLENGE SCHEMATIC .....                                                                                                                                        | 17        |
| 12.2      | DIFFUSIONAL FLOW CORRELATION .....                                                                                                                                                | 17        |
| 12.3      | WATER INTRUSION CORRELATION .....                                                                                                                                                 | 19        |
| <b>13</b> | <b>RETENTION TO AEROSOLISED <i>BREVUNDIMONAS DIMINUTA</i>.....</b>                                                                                                                | <b>20</b> |
| <b>14</b> | <b>CHEMICAL COMPATIBILITY .....</b>                                                                                                                                               | <b>21</b> |
| 14.1      | CHEMICAL COMPATIBILITY SUMMARY CHART FOR PHARMACEUTICAL PRODUCTS .....                                                                                                            | 21        |
| 14.2      | CHEMICAL COMPATIBILITY USER INSTRUCTIONS AND NOTES .....                                                                                                                          | 23        |

## 1 Introduction

Sterilisation grade filters used to sterilise gases that come into contact with pharmaceutical products must conform to strictly defined quality standards.

By using filter technology that conforms to the standards laid down by the various certifying bodies, the quality of the final product can be assured. Contamination can also be prevented from entering the final product by its comprehensive removal at each stage of the primary and secondary process.

When sterilising grade filters are used in the manufacture of products, the interactions between product, filter and process must be fully investigated and validated.

Guidelines for validation can be sourced from publications issued by the FDA, EMEA, USP, EP, BP, PDA<sup>1</sup> etc. This Validation Document has been produced with these guidelines in mind to enable the end user to incorporate this information within their own validation documentation or standard operating instructions for the process.

This Validation Document shows that TETPOR AIR surpasses the product specification requirements that have to be imposed on sterilising grade filters.

### NOTE

This Validation Document, if printed with the banner **UNCONTROLLED COPY**, has not been registered. If you wish to receive amendments and updates automatically, then please contact the department below or fax this page with your name and address and we will send you a registered copy: -

The Quality Assurance Department  
**domnick hunter ltd.**

#### **PROCESS FILTRATION**

Durham Road, Birtley  
Co. Durham, England, UK DH3 2SF  
Tel: +44 (0)191 4105121  
Fax: +44 (0)191 4105312  
E-mail: [process@domnickhunter.com](mailto:process@domnickhunter.com)

### **Registered Holder Information**

Name: .....

Company: .....

Department: .....

Address: .....  
.....  
.....

Registered No. ....

Date Issued: .....

Signed .....

<sup>1</sup> FDA, EMEA, USP, EP, BP, PDA – Food and Drug Administration, European Medicines Evaluation Agency, United States, European, and British Pharmacopoeia, Parenteral Drugs Association.

## 2 Quality Policy



### 2.1 Quality Assurance

Quality is of paramount importance to domnick hunter ltd. The filtration products are manufactured under controlled environmental conditions to the highest quality regimes and are subjected to a demanding program of Quality Assurance. Inspection and test protocols are implemented from vendor assessment, specification and receipt of raw materials, through every stage of the manufacturing process culminating in a non-destructive integrity test of the filter prior to packing and release.

Every stage of the manufacturing process has well defined assembly protocols laid down thus ensuring operational repeatability.

**domnick hunter's** responsibility as a manufacturer of quality extends beyond the site at Birtley in the UK, through to the world-wide network of filter and separation specialists.


The Quality Assurance Department operates with a well-equipped Laboratory Services Department, in which specialised personnel are employed to perform the essential quality inspections. **domnick hunter** has been assessed by the British Standards Institution and is registered to BS EN ISO 9001:1994, which defines the standards for quality systems, model for quality assurance in design, development, production, installation and servicing. As a concerned environmentally aware manufacturer, **domnick hunter** is also audited and certified to BS EN ISO 14001:1996 for its Environmental Management System.



### 3 Product Traceability

All pharmaceutical grade filters are non-destructively integrity tested and flushed with Purified Water. They are then dried using HEPA filtered air and sealed in a protective polyethylene bag within the controlled manufacturing environment prior to despatch.

To enable full traceability of all pharmaceutical grade filter products, each filter module is marked with an individual serial number. In addition, each filter product is marked with a lot number, product code and general description which is also shown on both the protective polyethylene bag in which the filter is sealed and on the outer surface of the final product packaging.



## 4 Summary of Test Methods

A number of the critical procedures used within this validation document are based on external methods sourced from publications issued by the USP, EP, FDA, ISO<sup>2</sup> and ASTM<sup>3</sup>. The remaining procedures are in-house methods prepared using **domnick hunter's** relevant experience and industry accepted standards in the fine filtration field. The comprehensive nature of these procedures ensures high batch to batch compliance to the specification.

| Membrane Validation |                        | Filter Validation and Testing |                            |                 |                                     |                        |
|---------------------|------------------------|-------------------------------|----------------------------|-----------------|-------------------------------------|------------------------|
| Physical            | Chemical / Biological  | Integrity Testing             | Retention                  | Flow Parameters | Durability                          | Chemical / Biological  |
| Bubble Point        | Chemical Compatibility | Bubble Point                  | Liquid Bacterial Challenge | Clean Air Flow  | Max Allowable Steam life            | Chemical Compatibility |
| Water / IPA Flow    | Extractables           | Pressure Hold                 | Particle / Fibre Shedding  |                 | Max continuous running Temperature  | Extractables           |
| Thickness           | Pyrogen Extractables   | Diffusional Flow              |                            |                 | Max Allowable Differential Pressure | Bio-Compatibility      |
| Bacterial Challenge | Oxidisables            | Aerosol Challenge             |                            |                 |                                     | Pyrogen Extractables   |
| Burst Strength      | Bio-Compatibility      | Water Intrusion               |                            |                 |                                     | Oxidisables            |
| Autoclavability     |                        |                               |                            |                 |                                     |                        |
| Pore Size Analysis  |                        |                               |                            |                 |                                     |                        |

<sup>2</sup> ISO – International Organisation for Standardisation

<sup>3</sup> ASTM – American Society for Testing and Materials

## 5 Preparation of Pharmaceutical Grade Filters

### 5.1 *Pharmaceutical Grade Standard*

TETPOR AIR filters must meet stringent standards to be certified P-grade product by **domnick hunter ltd.** The standards that must be met are:-

- Conformation to the requirements for non-fibre releasing filters as laid down in the United States Food and Drug Administration Regulations 21CFR211.72 and 210.3(b), (6).
- Effluent quality following a purified water flush must also be met as determined by the following tests:
  - ◇ Test for oxidisable material per USP 23 Purified Water.
  - ◇ Test for bacterial endotoxins using a gel clot LAL (Limulus Amoebocyte Lysate).
  - ◇ Test for particulates.
  - ◇ Test for TOC (Total Organic Carbon).
  - ◇ Test for Conductivity.
- All components conform to the Biological Safety Standards identified in USP <88> to Class VI-121°C levels.
- All filters must meet minimum defined integrity test values prior to despatch, defined from a correlation against an appropriate bacterial challenge.
- All filters are flushed with a high flux of purified water prior to despatch, as a guarantee of product cleanliness.

### 5.2 *Quality of Purified Water used in the preparation of Pharmaceutical Grade Filters*

The current USP and EP standards for Purified Water and Highly Purified Water specify a maximum conductivity of 1.1  $\mu\text{S}/\text{cm}$  @ 20°C (68°F) and a maximum TOC (Total Organic Carbon) content of 0.5 mg (500ppb) of carbon per litre.

The water used in the flushing stages of pharmaceutical grade filters exceeds these Pharmacopoeial requirements.

In addition, two other tests are carried out on samples taken from point of use and from a number of points in the supply pipework: -

- A Membrane Filtration Method Standard Plate Count Technique, to establish the microbial content.
- A gel clot LAL (Limulus Amoebocyte Lysate) test, for detection of bacterial endotoxins at 0.125 EU/ml sensitivity.

## 6 Certificate of Conformance

To certify that domnick hunter's TETPOR AIR filter products meet the highest pharmaceutical quality and performance requirements, a Certificate of Conformance is issued.

The example following is for TETPOR AIR cartridge filters.<sup>4</sup>


| Certificate of Conformance                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                              |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| For Pharmaceutical Grade                                                                                                                                                                                                                                                                                                                                         | Filters                                                                                                                                                                                      |
| Rated To                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                              |
| <i>This certifies that the domnick hunter filter</i>                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                              |
| <b>Recorded Lot Number</b>                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                              |
| has been manufactured in a purpose built facility within a controlled environment and subjected to a purified water flush.                                                                                                                                                                                                                                       |                                                                                                                                                                                              |
| <b>Materials of Construction</b>                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                              |
| All components of this cartridge are manufactured from materials suitable for contact with food and conform to the biological safety requirements as laid down in the current USP Class VI - 121°C Plastics. They also conform with the requirements for non fibre releasing filters as laid down in the United States FDA Title 21 CFR 211.72 and 210.3(b),(6). |                                                                                                                                                                                              |
| The filters also meet the domnick hunter quality control and assurance standards.                                                                                                                                                                                                                                                                                |                                                                                                                                                                                              |
| <b>Product Integrity</b>                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                              |
| This product has successfully passed a non destructive integrity test, which for sterilising grade cartridges is correlated to a bacterial challenge which gave a sterile effluent when challenged with a minimum of 10 <sup>7</sup> organisms per sq cm.                                                                                                        |                                                                                                                                                                                              |
| Diffusional flow value did not exceed      mL/min at a test pressure of      mbar                                                                                                                                                                                                                                                                                |                                                                                                                                                                                              |
| During validation, filter samples underwent the following tests and satisfactorily met the respective criteria specified:-                                                                                                                                                                                                                                       |                                                                                                                                                                                              |
| <b>Effluent Quality</b>                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                              |
| <b>TOC</b>                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                              |
| Met the requirements of USP Total Organic Carbon <643>                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                              |
| <b>Bacterial Endotoxins</b>                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                              |
| Cartridge aqueous extraction contains less than      EU/mL as determined using the Limulus Amebocyte Lysate (LAL) test, which meets the requirements of USP Water for Injection.                                                                                                                                                                                 |                                                                                                                                                                                              |
| <b>Water Conductivity</b>                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                              |
| Met the requirements of USP Water Conductivity <645>                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                              |
| <b>Particle Release</b>                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                              |
| Met the requirements of USP Particulate Matter in Injections <788>                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                              |
| <b>Thermal Stability</b>                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                              |
| Integrity was maintained after      steam cycles of      minutes at      °C                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                              |
| <br>FM 39187                                                                                                                                                                                                                                                                  | <p>This Product is registered with the Food &amp; Drug Administration<br/>Drug Master File No 7564</p>  |
| domnick hunter limited - Process Filtration, Durham Road, Birtley, Co. Durham, England DH3 2SF<br>Tel: +44 (0) 191 4105121 Fax: +44 (0) 191 4105312                                                                                                                                                                                                              |                                                                                                                                                                                              |

<sup>4</sup> A specific Certificate of Conformance is issued for DEMICAP capsules.

## 7 Product Coding and Range for TETPOR AIR Cartridges

Represented below are details of the product code structure for TETPOR AIR cartridges. This product code structure indicates cartridge sizes, micron ratings, endcap configurations and O'rings that are available within the product range.


Example: ZCMT2/-020C-A 500mm (20") 0.01 micron TETPOR AIR filter cartridge, pharmaceutical grade with 'C' style endcap and Silicone O'rings.



## 8 Product Coding and Range for TETPOR AIR Capsule Filters

Represented below are details of the product code structure for TETPOR AIR DEMICAP capsules. This product code structure indicates capsule sizes, micron ratings, connection configurations, sterility and material options that are available within the product range.

Example: ZEMTB-020TT-PN3 B-size 0.01 micron TETPOR AIR filter capsule, pharmaceutical grade with Tri-clamp connections, supplied non sterile in packs of 3.



<sup>5</sup> Silicone O'rings are fitted as standard without having to specify the S in the code.


<sup>6</sup> EPDM – Ethylene Propylene Diene Monomer Rubber

<sup>7</sup> Viton is a registered trademark of DuPont Dow Corporation


## 9 Endcap Configuration Diagrams




A Style  
223 'O' Rings




B,L Style  
Flat Gaskets



C Style  
226 'O' Rings




D Style  
222 'O' Rings



E Style  
222 'O' Rings



F Style  
216/218 'O' Rings



G Style  
222 'O' Rings




J Style  
S.O.E.



K Style  
214 'O' Rings




M Style  
214 'O' Rings




N Style  
213 'O' Rings


P Style  
227 'O' Rings



Q Style  
226 'O' Rings



R Style  
222 'O' Rings



SK Style



T Style  
126 'O' Rings



V Style  
Flat Gaskets



W Style  
111 'O' Rings



X Style  
Flat Gaskets



Y Style  
116 'O' Rings (Internal)




Z Style  
116 'O' Rings (Internal)




## 10 Capsule Dimensions


A Size




B Size




E Size





Fitted with 1 Tri-Clamp connections. Use drawings below to calculate dimensions with different configurations



1 Tri-Clamp



Stepped Hosebarb  
Suitable for tubing with 6mm ( $\frac{1}{4}$ ) -  
12mm ( $\frac{1}{2}$ ) internal diameter

 $\frac{1}{4}$  NPTM Thread $\frac{1}{2}$  Hosebarb $\frac{1}{2}$  NPTM ThreadFilling Bell with  
Stepped Hosebarb

NB. All dimensions nominal

## 11 Product Specification

### 11.1 Application

All products within the TETPOR AIR range, cartridge, capsule and disc have been designed for use in pharmaceutical applications for sterilisation duties, evaluated using ASTM Standard Test Method F838-83, or for prefiltration applications.

### 11.2 Materials of Construction for TETPOR AIR Cartridges and Capsules



All materials used in the construction of TETPOR AIR products that have product contact have met the requirements of the current USP Biological Reactivity Tests, In Vivo to Plastics Class VI-121°C.

All jointed surfaces are assembled by the use of heat sealing technology. No resins or binders are used in the manufacture of the filter and no surfactants are added to aid wetting.

| Component             | Material of Construction                     |
|-----------------------|----------------------------------------------|
| Filtration Membrane   | Polytetrafluoroethylene (PTFE)               |
| Upstream Support      | Polypropylene                                |
| Downstream Support    | Polypropylene                                |
| Inner Support Core    | Polypropylene                                |
| Outer Protection Cage | Polypropylene                                |
| Endcaps               | Polypropylene                                |
| Endcap Insert         | 316 Stainless Steel / Polysulphone (T-Style) |
| Standard O'rings      | Silicone                                     |
| Capsule Body          | Polypropylene                                |
| Capsule Vent Seals    | Silicone                                     |
| Filling Bell          | Polycarbonate                                |

### 11.3 Flow Characteristics

The effective filtration area of a standard 250mm (10") module is 0.75 m<sup>2</sup> (7.9 ft<sup>2</sup>).



### 11.4 Operating Temperatures and Pressures

Below are the recommended maximum differential operating pressures at various temperatures: -

| Temperature |     | Differential Pressure Cartridges |      | Operating Pressure Capsules<br>Liquids (gases) |             |
|-------------|-----|----------------------------------|------|------------------------------------------------|-------------|
| °C          | °F  | bar                              | psi  | barg                                           | psig        |
| 20          | 68  | 5.00                             | 72.5 | 5.00 (4.00)                                    | 72.5 (58.0) |
| 40          | 104 | 4.00                             | 58.0 | 5.00 (4.00)                                    | 72.5 (58.0) |
| 50          | 122 | 3.00                             | 43.5 | 3.00 (3.00)                                    | 43.5 (43.5) |
| 70          | 158 | 2.00                             | 29.0 | Not Recommended                                |             |

### 11.5 Steam Sterilisation

#### Autoclave

| Product Format | Autoclave Temp |     | Number of Cycles | Cycle Time (minutes at temperature) |
|----------------|----------------|-----|------------------|-------------------------------------|
|                | °C             | °F  |                  |                                     |
| Cartridges     | 142            | 288 | 120              | 30                                  |
| Capsules       | 135            | 275 | 10               | 30                                  |

To maximise cartridge and capsule life, a slow exhaust cycle is recommended.

#### Steam in Place (SIP)

| Product Format | SIP Temp |     | Number of Cycles | Cycle Time (minutes at temperature) |
|----------------|----------|-----|------------------|-------------------------------------|
|                | °C       | °F  |                  |                                     |
| Cartridges     | 142      | 288 | 120              | 30                                  |
| Capsules       | 135      | 275 | 10               | 30                                  |

To maximise the life of the cartridge, the differential pressure across the cartridge should not exceed 0.30 bar (4.4 psi) at 142°C (288°F). For new applications it is recommended that the **domnick hunter** guidance for the method of steam sterilisation be followed.

### 11.6 Tests For BioCompatibility

An independent research establishment has assessed the biological safety associated with the use of TETPOR AIR filters designed for processing pharmaceutical products. The materials used in the construction of TETPOR AIR products meet the requirements of the current USP <88> *Biological Reactivity* tests at Plastics Class VI – 121°C.

### 11.7 Extractables

All pharmaceutical grade filters are designed and manufactured to yield a minimum of extractables. Testing of a purified water filtrate with TETPOR AIR is documented below.

#### Test Method (1)

Non-volatile extractables from purified water samples after flowing through an autoclaved 250mm (10") TETPOR AIR cartridge are listed below. The levels shown are the quantities present in 100ml samples of filtrate which were taken at stages throughout a 10-litre flush.

**Cartridge Serial No. S000136**

| Flush Quantity (litres) | Non-volatile Extract (mg per 100ml) | Oxidisables (per USP 23 Test Method) |
|-------------------------|-------------------------------------|--------------------------------------|
| 1                       | 0.5                                 | PASS                                 |
| 2                       | 0.1                                 | PASS                                 |
| 4                       | <0.1                                | PASS                                 |
| 6                       | 0.1                                 | PASS                                 |
| 8                       | <0.1                                | PASS                                 |
| 10                      | <0.1                                | PASS                                 |

#### Test Method (2)

Non-volatile extractables from an autoclaved 250mm (10") TETPOR AIR cartridge and a B-size TETPOR AIR capsule following a 4-hour dynamic immersion in a variety of commonly used solvents. Solvent volume used 1500ml.

| Solvent                   | Serial No. Cartridge | Weight of extract (mg) | Serial No. Capsule | Weight of extract (mg) |
|---------------------------|----------------------|------------------------|--------------------|------------------------|
| Water                     | M41711               | 2.2                    | A022783            | <0.01                  |
| Iso Propyl Alcohol (IPA)  | M41719               | 17.6                   | A022782            | 0.38                   |
| Ethanol                   | M41729               | 9.8                    | A022777            | 0.11                   |
| Methanol                  | M41717               | 12.0                   | A022773            | 3.00                   |
| Ammonium Hydroxide (pH12) | M41708               | 0.8                    | A022776            | <0.01                  |
| Hydrochloric Acid (pH3)   | M41714               | 8.3                    | A022774            | <0.01                  |

## 11.8 Integrity Testing Data

The following is the integrity test information for the micron ratings available within the TETPOR AIR product range. Diffusional flow and bubble point values are given for cartridges wetted in 60:40<sub>v/v</sub> IPA:Water solution.

| Micron Rating | Minimum Bubble Point <sup>8</sup> |      | Diffusional Flow Test Pressure |      | Maximum Diffusional Flow (ml/min)  |     |                |                |                |
|---------------|-----------------------------------|------|--------------------------------|------|------------------------------------|-----|----------------|----------------|----------------|
|               | bar                               | psi  | bar                            | psi  | 10"                                | K   | A <sup>9</sup> | B <sup>7</sup> | E <sup>7</sup> |
| Liquid / Gas  |                                   |      |                                |      |                                    |     |                |                |                |
| 0.20 / 0.01   | 1.00                              | 14.5 | 0.80                           | 11.6 | 18.0                               | 8.5 | 6.0            | 3.0            | 1.5            |
|               |                                   |      | Water Intrusion Test Pressure  |      | Maximum Water Intrusion (ml/10min) |     |                |                |                |
|               |                                   |      | bar                            | psi  | 10"                                | K   | A <sup>7</sup> | B <sup>7</sup> | E <sup>7</sup> |
|               |                                   |      | 2.50                           | 36.3 | 16.0                               | 7.5 | 5.3            | 2.6            | 1.3            |

## 12 Correlation of Non-Destructive Integrity Testing to Liquid Bacterial Challenge with *Brevundimonas diminuta* (ATCC 19146) for 0.20 µm Sterilising Grade Filters

To guarantee filter performance a filter must be capable of being non-destructively integrity tested. This was recognised by the FDA in the "Guideline on Sterile Drug Products produced by Aseptic Processing" (June 1987), which states:

*"After a filtration process is properly validated for a given product process and filter, it is important to assure that identical filter replacements (membrane or cartridge) used in production runs will perform in the same manner. One way of achieving this is to correlate filter performance data with filter integrity testing data. Normally, integrity testing of the filter is performed after the filter unit is assembled and sterilised prior to use. More importantly, however, such testing should be conducted after the filter is used in order to detect any filter leaks or perforations that may have occurred during the filtration."*


To achieve this objective the correlation between bacterial challenge retention and a non-destructive integrity test must be proven. The procedure documented in ASTM F838-83, 'Standard Test Method for Determining Bacterial Retention Of Membrane Filters Utilized For Liquid Filtration', is used to test the manufactured product<sup>10</sup>. The filter must be challenged with a minimum of 10<sup>7</sup> viable *Brevundimonas diminuta* (ATCC 19146) per cm<sup>2</sup> of effective filtration area.

<sup>8</sup> domnick hunter does not recommend the use of bubble point as an integrity test method for cartridges and capsules, but values are given for use as an indicator of product integrity.

<sup>9</sup> DEMI Cartridge and DEMICAP Capsule.

<sup>10</sup> Previous reference to the guidance document *Microbial Evaluation of Filters for Sterilising Liquids*, HIMA Document No. 3 Vol. 4, April 1982, referred to in USP<1211> *Sterilisation by Filtration* has been superseded by the equivalent ASTM F838-83.

## 12.1 Liquid Bacterial Challenge Schematic



Any organisms that pass through the test filter are collected and cultured on the surface of analytical discs. In this way colonies may be counted and bacterial species identified. The bacterial challenge is quantified by expressing the filters' efficiency to remove the challenge organism from the challenge suspension as a Log Reduction Value (LRV). To be classed as a sterilising grade filter the test filter must produce a sterile filtrate and  $LRV > 7.0$ .

$$LRV = \log_{10} \left( \frac{\text{Number of organisms in the challenge}}{\text{Number of organisms in the filtrate}} \right)$$

## 12.2 Diffusional Flow Correlation

Table 1 lists TETPOR AIR 0.01  $\mu\text{m}$  cartridges that were diffusional flow integrity tested before and after bacterial challenge. The bacterial challenge was conducted using ASTM F838-83 so providing the necessary correlation between a bacterial challenge and a non-destructive diffusional flow test.

Table 1 indicates that a 250mm (10") TETPOR AIR 0.01  $\mu\text{m}$  filter exhibiting a diffusional flow of  $<22.0 \text{ ml/min}$  when completely wet out in 60:40<sub>v/v</sub> IPA:Water, at a test pressure of 0.80 barg (11.6 psig) at 20°C (68°F) will produce a sterile filtrate. To build in an added level of security, **domnick hunter** has reduced the allowable diffusional flow.

**A 250mm (10") TETPOR AIR 0.01  $\mu\text{m}$  filter, that is completely wet out in 60:40<sub>v/v</sub> IPA:Water, has a maximum diffusional flow value of:**

**18.0 ml/min at a test pressure of 0.80 barg (11.6 psig) at 20°C (68°F).**

The maximum diffusional flow for other TETPOR AIR filters is based upon the ratio of the surface area compared to that of a 250mm (10") cartridge.

**Table 1: Diffusional Flow Correlation Data**

Filter Type: ZCMT1/-020C-A TETPOR AIR 0.01µm 10" Cartridge

Challenge Organism: *Brevundimonas diminuta* (ATCC 19146)

| Serial No. | Diffusional Flow<br>(Air in 60:40 <sub>v/v</sub><br>IPA:Water) @ 0.80<br>barg (11.6 psig)<br>ml/min | Challenge<br>Level (x 10 <sup>11</sup> ) | Organisms<br>Passed | LRV <sup>11</sup> |
|------------|-----------------------------------------------------------------------------------------------------|------------------------------------------|---------------------|-------------------|
| M32898     | 0.5                                                                                                 | 1.70                                     | 0                   | 11.23             |
| M36067     | 2.0                                                                                                 | 2.50                                     | 0                   | 11.40             |
| M32903     | 2.2                                                                                                 | 1.34                                     | 0                   | 10.64             |
| M32914     | 2.2                                                                                                 | 2.64                                     | 0                   | 10.21             |
| M36065     | 2.3                                                                                                 | 0.31                                     | 0                   | 10.49             |
| M32894     | 3.1                                                                                                 | 1.64                                     | 0                   | 11.21             |
| M32902     | 3.1                                                                                                 | 1.77                                     | 0                   | 11.25             |
| M41707     | 3.5                                                                                                 | 0.14                                     | 0                   | 10.10             |
| M32908     | 3.7                                                                                                 | 2.41                                     | 0                   | 11.38             |
| M41710     | 5.1                                                                                                 | 0.15                                     | 0                   | 10.17             |
| M36069     | 6.1                                                                                                 | 1.57                                     | 0                   | 11.20             |
| N008205    | 6.6                                                                                                 | 2.88                                     | 0                   | 11.45             |
| M37703     | 6.7                                                                                                 | 3.25                                     | 0                   | 11.51             |
| M32916     | 6.7                                                                                                 | 2.32                                     | 0                   | 11.36             |
| M32910     | 6.8                                                                                                 | 1.65                                     | 0                   | 11.22             |
| M41722     | 6.8                                                                                                 | 0.13                                     | 0                   | 10.11             |
| N008204    | 7.0                                                                                                 | 1.99                                     | 0                   | 11.30             |
| M37700     | 7.2                                                                                                 | 3.08                                     | 0                   | 11.40             |
| J29578     | 7.9                                                                                                 | 1.80                                     | 0                   | 11.18             |
| J29577     | 7.9                                                                                                 | 3.02                                     | 0                   | 11.48             |
| S000143    | 9.9                                                                                                 | 2.87                                     | 0                   | 11.46             |
| M36064     | 11.0                                                                                                | 3.11                                     | 0                   | 11.49             |
| M37710     | 11.0                                                                                                | 2.01                                     | 0                   | 11.30             |
| E92447     | 22.0                                                                                                | 0.52                                     | 0                   | 10.72             |
| M38639     | 22.4                                                                                                | 3.29                                     | 15                  | 10.30             |
| M41723     | 36.5                                                                                                | 0.40                                     | 333                 | 7.79              |

<sup>11</sup> Where Organisms passed = 0, LRV is stated as *greater than*.

### 12.3 Water Intrusion Correlation

Table 2 proves the correlation that a 250mm (10") TETPOR AIR 0.01um filter with a Water Intrusion value of <20.0 ml / 10min at a test pressure of 2.50 barg (36.3 psig) at 20°C (68°F) will produce a sterile filtrate. To build in an added level of security, **domnick hunter** has reduced the allowable Water Intrusion value.

**A 250mm (10") TETPOR AIR 0.01 um filter has a maximum allowable Water Intrusion value of:**

**16.0 ml / 10min at a test pressure of 2.50 barg (36.3 psig) at 20°C (68°F).**

The Water Intrusion values for other TETPOR AIR filters are based upon the ratio of the surface area compared to that of a 250mm (10") cartridge.

**Table 2: Water Intrusion Correlation Data**

Filter Type: ZCMT1/020C-A TETPOR AIR 0.01µm 10" Cartridge

Challenge Organism: *Brevundimonas diminuta* (ATCC 19146)

| Serial No. | WI @ 2.50 barg<br>(36.3 psig)<br>ml/10min | Challenge<br>Level (x 10 <sup>11</sup> ) | Organisms<br>Passed | LRV <sup>12</sup> |
|------------|-------------------------------------------|------------------------------------------|---------------------|-------------------|
| M36064     | 6.0                                       | 3.11                                     | 0                   | 11.49             |
| M36067     | 7.7                                       | 2.50                                     | 0                   | 11.39             |
| M36069     | 7.9                                       | 1.57                                     | 0                   | 11.19             |
| M32894     | 10.3                                      | 1.64                                     | 0                   | 11.21             |
| M32898     | 11.2                                      | 1.70                                     | 0                   | 11.23             |
| M32907     | 11.8                                      | 2.04                                     | 0                   | 11.31             |
| M32908     | 12.1                                      | 2.41                                     | 0                   | 11.38             |
| M32910     | 12.7                                      | 1.65                                     | 0                   | 11.22             |
| M32902     | 13.3                                      | 1.77                                     | 0                   | 11.25             |
| N008204    | 13.3                                      | 1.99                                     | 0                   | 11.30             |
| M32896     | 15.4                                      | 1.93                                     | 0                   | 11.29             |
| M32916     | 16.0                                      | 2.32                                     | 0                   | 11.37             |
| M32915     | 20.1                                      | 1.99                                     | 0                   | 11.30             |
| N008200    | 29.3                                      | 1.19                                     | TNTC <sup>13</sup>  | <7.00             |
| M38638     | 38.4                                      | 1.62                                     | 2083                | 7.90              |
| M38636     | 74.2                                      | 3.29                                     | 16                  | 10.30             |

<sup>12</sup> Where Organisms passed = 0, LRV is stated as *greater than*.

<sup>13</sup> TNTC – Too Numerous To Count

### **13 Retention to Aerosolised *Brevundimonas Diminuta***

Tests have shown that TETPOR AIR is fully retentive to aerosolised *Brevundimonas diminuta* (ATCC 19146) bacteria when challenged with a total of  $2 \times 10^9$  cfu over a 1-hr test at the rated flow of the cartridge.

## 14 Chemical Compatibility

Testing has been carried out that shows TETPOR AIR filters and capsules have a broad range of chemical compatibility with chemicals commonly used in the pharmaceutical industry.

### 14.1 Chemical Compatibility Summary Chart for Pharmaceutical Products

|                                | ASYPOR | BIO-X II | HIGH FLOW BIO-X | HIGH FLOW BIO-X VENT AUTOCLAVE | HIGH FLOW PREPOR GFA | HIGH FLOW TETPOR | HIGH FLOW TETPOR H.T. | HIGH FLOW TETPOR VENT AUTOCLAVE | PEPLYN PLUS & PEPLYN PLUS DC | PREPOR GF | PREPOR PES | PREPOR PES / FS | TETPOR AIR | TETPOR LIQUID | TETPOR PLUS | EPDM O'Ring | VITON O'Ring | SILICONE O'Ring |
|--------------------------------|--------|----------|-----------------|--------------------------------|----------------------|------------------|-----------------------|---------------------------------|------------------------------|-----------|------------|-----------------|------------|---------------|-------------|-------------|--------------|-----------------|
| Acetic acid 3.5N               | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Acetic acid 8.75N              | LC     | C        | C               | -                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | LC          | LC           | NC              |
| Acetic acid conc.17.5N         | NC     | C        | C               | -                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | LC          | NC           | NC              |
| Acetone                        | NC     | C        | C               | -                              | C                    | C                | C                     | C                               | C                            | C         | NC         | NC              | C          | C             | C           | NC          | NC           | NC              |
| Acetonitrile                   | NC     | C        | C               | -                              | LC                   | C                | C                     | C                               | C                            | LC        | -          | -               | C          | C             | C           | NC          | NC           | NC              |
| Acidbrite 4 (Diversey) 3.0%v/v | NC     | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | -          | -               | -          | -             | -           | C           | C            | C               |
| Ammonium Hydroxide 8N          | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | LC              | C          | C             | C           | C           | C            | C               |
| Ammonium Oxalate 0.07N         | -      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Amyl Acetate                   | NC     | C        | C               | C                              | LC                   | C                | C                     | C                               | C                            | LC        | LC         | LC              | C          | C             | C           | NC          | NC           | LC              |
| Aqueous Ammonia 15.5N          | NC     | C        | C               | C                              | LC                   | C                | LC                    | C                               | C                            | LC        | C          | LC              | C          | C             | C           | C           | C            | C               |
| Benzyl Alcohol                 | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | NC                           | NC        | -          | -               | C          | C             | C           | C           | C            | C               |
| Benzalkonium Chloride 0.1%     | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Boric acid,saturated           | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Butan-1-ol                     | NC     | C        | C               | C                              | C                    | LC               | LC                    | LC                              | C                            | C         | C          | C               | NC         | NC            | NC          | C           | C            | C               |
| Butan-2-ol                     | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | LC          | C            | C               |
| Carbon Tetrachloride           | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | NC                           | NC        | -          | -               | NC         | NC            | NC          | NC          | C            | NC              |
| Chloroform                     | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | NC                           | NC        | NC         | NC              | NC         | NC            | NC          | NC          | LC           | NC              |
| Cyclohexane                    | NC     | C        | C               | C                              | NC                   | -                | -                     | -                               | NC                           | NC        | -          | -               | LC         | LC            | LC          | NC          | NC           | NC              |
| 1,4 – Dioxane                  | NC     | C        | C               | C                              | LC                   | C                | C                     | C                               | C                            | LC        | -          | -               | C          | C             | C           | NC          | NC           | NC              |
| Diverflow (Diversey) 3%v/v     | NC     | -        | -               | -                              | NC                   | -                | -                     | -                               | C                            | NC        | C          | C               | -          | -             | -           | C           | C            | LC              |
| Diversey 212G 0.6%v/v          | NC     | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | -          | -               | -          | -             | -           | C           | C            | C               |
| Divosan Forte 0.5%v/v          | LC     | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | C          | C               | -          | -             | -           | C           | C            | C               |
| Divosan XT 1%v/v               | C      | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | -          | -               | -          | -             | -           | C           | C            | C               |
| Ethanol                        | NC     | C        | C               | C                              | C                    | C                | -                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | LC              |
| Ethanol 45%                    | LC     | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Ethyl Acetate                  | NC     | LC       | LC              | LC                             | LC                   | LC               | LC                    | LC                              | LC                           | LC        | NC         | NC              | LC         | LC            | LC          | C           | NC           | LC              |
| Formaldehyde 0.3%              | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Formaldehyde 37%               | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Formic acid conc.              | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | -          | -               | C          | C             | C           | C           | NC           | NC              |

|                                     | ASYPOR | BIO-X II | HIGH FLOW BIO-X | HIGH FLOW BIO-X VENT AUTOCLAVE | HIGH FLOW PREPOR GFA | HIGH FLOW TETPOR | HIGH FLOW TETPOR H.T. | HIGH FLOW TETPOR VENT AUTOCLAVE | PEPLYN PLUS & PEPLYN PLUS DC | PREPOR GF | PREPOR PES | PREPOR PES / FS | TETPOR AIR | TETPOR LIQUID | TETPOR PLUS | EPDM O'Ring | VITON O'Ring | SILICONE O'Ring |
|-------------------------------------|--------|----------|-----------------|--------------------------------|----------------------|------------------|-----------------------|---------------------------------|------------------------------|-----------|------------|-----------------|------------|---------------|-------------|-------------|--------------|-----------------|
| Glycerol                            | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Hexane                              | LC     | C        | C               | C                              | -                    | C                | C                     | C                               | NC                           | -         | -          | -               | -          | -             | -           | NC          | NC           | NC              |
| Hydrochloric acid 1N                | C      | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Hydrochloric acid conc.             | NC     | -        | -               | -                              | NC                   | -                | -                     | -                               | C                            | NC        | -          | -               | C          | C             | C           | NC          | NC           | NC              |
| Hydrochloric acid conc.13%          | -      | C        | C               | C                              | -                    | C                | C                     | C                               | -                            | -         | -          | -               | -          | -             | -           | NC          | NC           | NC              |
| Hydrogen Peroxide                   | -      | C        | C               | C                              | -                    | -                | -                     | -                               | C                            | -         | -          | -               | -          | -             | -           | C           | C            | C               |
| Hydrogen Peroxide 10 Vol            | C      | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Hydrogen Peroxide 100 Vol           | LC     | -        | -               | -                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Methanol                            | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | NC           | C               |
| Methyl-Iso-Butylketone              | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | NC         | NC              | C          | C             | C           | NC          | NC           | LC              |
| Methylene Chloride @ 40°C           | -      | -        | -               | -                              | LC                   | -                | -                     | -                               | LC                           | LC        | -          | -               | -          | -             | -           | -           | -            | -               |
| Nitric acid 2N 14.4%                | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | LC          | C            | C               |
| Nitric acid 15.8N                   | NC     | C        | C               | C                              | NC                   | C                | NC                    | C                               | C                            | NC        | -          | -               | C          | C             | C           | NC          | NC           | NC              |
| Ozone                               | -      | -        | -               | -                              | -                    | -                | -                     | -                               | -                            | -         | NC         | NC              | -          | -             | -           | -           | -            | -               |
| Paraffin yellow                     | LC     | LC       | LC              | LC                             | LC                   | C                | C                     | C                               | C                            | LC        | -          | -               | C          | C             | C           | NC          | C            | NC              |
| Pentane                             | LC     | C        | C               | C                              | LC                   | -                | -                     | -                               | LC                           | LC        | -          | -               | LC         | LC            | LC          | NC          | C            | NC              |
| Peracetic acid 0.5% (10 wk test)    | C      | -        | -               | -                              | -                    | C                | C                     | C                               | -                            | -         | -          | -               | -          | -             | -           | C           | C            | C               |
| Peracetic acid 4%                   | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Perchloroethylene                   | -      | -        | -               | -                              | -                    | -                | -                     | -                               | -                            | -         | NC         | NC              | -          | -             | -           | -           | -            | -               |
| Petroleum spirits                   | NC     | -        | -               | -                              | NC                   | C                | C                     | C                               | NC                           | NC        | -          | -               | LC         | LC            | LC          | NC          | C            | NC              |
| Phenol (aq) 0.5N                    | -      | C        | C               | C                              | -                    | NC               | -                     | NC                              | -                            | -         | -          | -               | -          | -             | -           | -           | -            | -               |
| Phenol 5%                           | NC     | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Phenol 0.25%                        | C      | -        | -               | -                              | C                    | -                | -                     | -                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Polyethylene Glycol 600             | NC     | LC       | LC              | LC                             | NC                   | C                | C                     | C                               | LC                           | NC        | NC         | NC              | -          | -             | -           | -           | -            | -               |
| Polyglycol 2000-E                   | -      | -        | -               | -                              | -                    | -                | -                     | -                               | -                            | -         | -          | -               | -          | -             | -           | C           | C            | C               |
| Potassium Dichromate 0.1N           | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Potassium Iodine 0.6N               | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Potassium Hydroxide 10N             | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | C          | LC              | C          | C             | C           | C           | C            | C               |
| Potassium Permanganate 0.1N         | NC     | C        | C               | C                              | NC                   | C                | LC                    | C                               | C                            | NC        | C          | C               | C          | C             | C           | C           | C            | C               |
| Propan-1-ol                         | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | C          | C               | C          | C             | C           | C           | C            | LC              |
| Propan-2-ol                         | C      | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | C          | C               | C          | C             | C           | C           | C            | LC              |
| Propan-2-ol, 60:40 H <sub>2</sub> O | C      | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | C          | C               | C          | C             | C           | C           | C            | C               |
| Pyridine                            | NC     | C        | C               | C                              | NC                   | C                | C                     | C                               | C                            | NC        | NC         | NC              | C          | C             | C           | C           | NC           | C               |
| Sodium Chloride 0.5N                | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Saline Lactose Broth                | C      | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |

|                                                 | ASYPOR | BIO-X II | HIGH FLOW BIO-X | HIGH FLOW BIO-X VENT AUTOCLAVE | HIGH FLOW PREPOR GFA | HIGH FLOW TETPOR | HIGH FLOW TETPOR H.T. | HIGH FLOW TETPOR VENT AUTOCLAVE | PEPLYN PLUS & PEPLYN PLUS DC | PREPOR GF | PREPOR PES | PREPOR PES / FS | TETPOR AIR | TETPOR LIQUID | TETPOR PLUS | EPDM O'Ring | VITON O'Ring | SILICONE O'Ring |
|-------------------------------------------------|--------|----------|-----------------|--------------------------------|----------------------|------------------|-----------------------|---------------------------------|------------------------------|-----------|------------|-----------------|------------|---------------|-------------|-------------|--------------|-----------------|
| Sodium Hydroxide 2N 8%                          | NC     | NC       | NC              | NC                             | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Sodium Hydroxide 7N 28%                         | NC     | NC       | NC              | NC                             | NC                   | C                | C                     | C                               | C                            | NC        | NC         | NC              | C          | C             | C           | C           | C            | LC              |
| Sodium Hypochlorite (14% Free Cl <sub>2</sub> ) | NC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | C          | C               | C          | C             | C           | C           | C            | C               |
| Sodium thiosulphate 0.1N                        | LC     | C        | C               | C                              | C                    | C                | C                     | C                               | C                            | C         | -          | -               | C          | C             | C           | C           | C            | C               |
| Sulphuric acid 1N                               | NC     | C        | C               | C                              | LC                   | C                | C                     | C                               | C                            | LC        | C          | C               | -          | -             | -           | C           | C            | C               |
| Sulphuric acid conc.                            | NC     | NC       | NC              | NC                             | LC                   | LC               | NC                    | LC                              | LC                           | LC        | NC         | NC              | LC         | LC            | LC          | -           | -            | -               |
| Sulphurous acid                                 | -      | -        | -               | -                              | -                    | -                | -                     | -                               | -                            | -         | NC         | NC              | -          | -             | -           | -           | -            | -               |
| Toluene                                         | -      | NC       | NC              | NC                             | -                    | NC               | NC                    | NC                              | NC                           | -         | NC         | NC              | -          | -             | -           | NC          | LC           | NC              |
| 1,1,1 Trichloroethane                           | LC     | -        | -               | -                              | -                    | -                | -                     | -                               | -                            | -         | -          | -               | -          | -             | -           | -           | -            | -               |
| 1,1,2 Trichloroethane                           | -      | C        | C               | C                              | LC                   | C                | LC                    | C                               | LC                           | LC        | NC         | NC              | LC         | LC            | LC          | NC          | LC           | LC              |
| Trichloroacetic Acid 80%                        | LC     | -        | -               | -                              | LC                   | -                | -                     | -                               | C                            | LC        | -          | -               | C          | C             | C           | NC          | LC           | NC              |
| Trichloroacetic Acid 5N                         | -      | C        | C               | C                              | -                    | C                | C                     | C                               | -                            | -         | -          | -               | -          | -             | -           | ---         |              |                 |
| Toluene                                         | NC     | -        | -               | -                              | NC                   | -                | -                     | -                               | -                            | NC        | -          | -               | -          | -             | -           | NC          | LC           | NC              |
| Xylene                                          | NC     | LC       | LC              | LC                             | NC                   | LC               | LC                    | LC                              | NC                           | NC        | LC         | LC              | NC         | NC            | NC          | C           | LC           | NC              |

## 14.2 Chemical Compatibility User Instructions and Notes

- The chemicals are arranged in alphabetical order using their most common or trade names. If the chemical in question does not appear to be listed, it may be found elsewhere in the table under a pseudonym, in particular its IUPAC<sup>14</sup> name.
- With regard to compatibility:
  - ◇ Any product that has Limited Compatibility (LC) at ambient temperatures should not be used at a higher temperature.
  - ◇ The list of compatibilities does not take into account any synergistic effects of more than one chemical present in the solution to be filtered.
  - ◇ Test Conditions – 72 hrs at ambient temperature and pressure, unless otherwise stated.
  - ◇ Contact **domnick hunter** for confirmation of compatibility with specific operating conditions.

<sup>14</sup> IUPAC – International Union of Pure and Applied Chemistry

**Documentation Approval Section**

Approved By: \_\_\_\_\_ Dr M Brailsford \_\_\_\_\_  
Title: \_\_\_\_\_ Technical Director \_\_\_\_\_  
Signature: \_\_\_\_\_  
Date: \_\_\_\_\_ 20 May 2003 \_\_\_\_\_



**domnick hunter**  
**PROCESS FILTRATION**

Durham Road, Birtley, Co. Durham  
England, UK. DH3 2SF  
Tel: +44 (0)191 4105121  
Fax: +44 (0)191 4105312  
E-mail: [process@domnickhunter.com](mailto:process@domnickhunter.com)  
WebSite: [www.domnickhunter.com](http://www.domnickhunter.com)

**domnick hunter limited** has a policy of continuous product development and although the Company reserves the right to change specification, it attempts to keep customers informed of any alterations. This publication provides generic validation information only and customers are requested to contact the Process Division for detailed information and advice on a product's suitability for specific applications.

Copyright **domnick hunter limited** 2003